

Document Title: Description LIT, SDS, Indicating Drierite		Part # and Rev. 18025	
		Release Date: 10/09/15	
Rev.	Description	Revision By	Date
1	Updated to new SDS from manufacturer. Clarified that this is indicating Drierite	JCP	4/18/19

Production Filename: 18025

Working Path: METER Product Number Library/18000-18099/18025

Dimensions: 5.5 inches wide, 8.5 inches tall

Material: Plain White Paper

Colors: Black printing

Printer Type: Any Laser Printer that can duplex

Special Notes: Fold in half

Illustrations are Ref Only ** Not to Scale **

W.A. Hammond Drierite Co., LTD

Safety Data Sheet

1. PRODUCT AND COMPANY IDENTIFICATION

1.1 Product Identifiers	
Product Name	Drierite, Indicating
Product Stock Numbers	21001, 21005, 21025, 22001, 22005, 22025, 22001, 23005, 23025, 24001, 24005, 24025, 24035, 26800, 26840, 26930, 27068, 27069, 27070, 40207, 40451, 50068
Manufacturer	W.A. Hammond Drierite Co., LTD.
1.2 Relevant identified uses of the substance or mixture and uses advised against	
Identified Uses	Desiccant or Drying Agent
1.3 Details of the supplier of the safety data sheet	
Company	W.A. Hammond Drierite Co., LTD. P.O. Box 460 Xenia, OH 45385 U.S.A.
Telephone	937-376-2927
Website	www.drierite.com
1.4 Emergency telephone number	
Emergency Phone#	937-376-2927

2. HAZARDS IDENTIFICATION

2.1 Classification of the substance or mixture	
GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)	
Eye Irritation (Category 2A), H319	
Respiratory Sensitization (Category 1), H334	
Skin Sensitization (Category 1), H317	
Germ Cell Mutagenicity (Category 2), H341	
Carcinogenicity (Category 2), H351	
Reproductive Toxicity (Category 1B), H360	
Acute Aquatic Toxicity (Category 2), H401	
Chronic Aquatic Toxicity (Category 3), H412	
For the full text of the H-Statements mentioned in this section, see Section 16.	
2.2 GHS Label elements, including precautionary statements	
Pictogram	
Signal Word	Danger
Hazard Statement(s)	H317
	May cause an allergic skin reaction.