

Document Title: Description Box, PTT Standards (All types)		Part # 13997 13998 13999 14500 14501 14502 18053	
		Release Date: 9/28/2011	
Rev.	Description	Revision By	Date
1	File Created	TJZ	9/28/11
	Uploaded sticker to repository in the Inserts folder. Available at http://publications.decagon.com . Please ask archivist for previous versions or use Beanstalk application.	DDH	2/7/14

Production File: [http://manuals.decagon.com/Inserts/ Verification Standards Box Art_Print.eps](http://manuals.decagon.com/Inserts/VerificationStandardsBoxArt_Print.eps)

Die Line file (Ai & PDF)

http://manuals.decagon.com/Inserts/UNICEP_DISPENSER_BOX - 4-21-2011.ai

http://manuals.decagon.com/Inserts/UNICEP_DISPENSER_BOX - 4-21-2011.pdf

Size 3.5x2.75x4.25

Materials .020 SBS

Graphics 4 Colors + 1 spot color (0.25 Red Part #13997, 0.50 Green Part #13998, 0.76 Blue Part #13999, .920 Purple Part #14500, 0.98 Orange Part #14501, 1.000 Lite Blue Part #14502, .150 Cool Gray 10C #18053.)

Perforation: Pull Tab

Special Notes:

Please make sure that the box folds together square and the pull tab works fine.

All types of standards box are in a single file with separate layers. The first layer is the die line, second layer is the CMYK plate and the other layers correspond with their individual part numbers.

Illustrations are for Ref Only ****Not to Scale****

PRINT VIEW

CAD #5941

0.760 a_w
NaCl 6.00 mol/kg in H₂O

Temperature Correction of Decagon's Verification Standards

Temp °C	6.00 mol/kg NaCl a_w
15.0	0.761
20.0	0.760
25.0	0.760
30.0	0.760
35.0	0.760
40.0	0.760

Reorder at www.aqualab.com/reorder
For orders outside of the USA, contact your local distributor.
Phone: 1-800-755-2751
email: support@aqualab.com



Time for a factory calibration?

- Accurate readings are important to me.
- I want independent or factory verification.
- My calibration certificate has expired.
- I want to extend my warranty.


Decagon recommends preventative maintenance plans and annual calibration service for your AquaLab. Call or visit our website to learn how automatic service plans make it easy.

AQUA LAB
BY DECAGON

COUNT: 50 VIALS

- ✓ Machine verified standard within specifications
- ✗ Machine did not verify standard within specifications

Warm the machine up for at least 15 minutes.



*See the verification section in the operator's manual

Visit www.aqualab.com/training for a video demonstration and other tips.