

Decagon Devices, Inc. Product Market Information Sheet

Product: Decagon ProCheck

Description: Instantaneous readout device for all of Decagon's sensors. The ProCheck can store 5,000 readings including sensor type, date, time, raw value, calibrated values, and calibration coefficients. The ProCheck connects to computers using Serial RS232 for easy data download. The ProCheck is powered by 4 AA batteries.

Distributors of Product: Decagon Devices, Inc. is the sole distributor of the ProCheck in the USA to universities and other research institutions.

Competitors with similar products: None known with the capabilities of the ProCheck.